

## BRG Guideline: Basic NMR Spectrometer Operation

The following guidelines provide an overview to operation of the departmental Bruker NMR spectrometers and are to be used as a memory aid for new users. Before using NMR spectrometers for the first time, users **must** attend a training session. Rodger Kohnert should be consulted immediately if any problems are encountered when using the instruments.

### 1. Arrival at NMR

- log-in and run Topspin software (single click on icon)
- press "lift" button to engage sample air lift
- place NMR sample tube in spinner and adjust to correct depth with gauge
- gently place NMR tube/spinner on air cushion in center of magnet
- turn-off air lift (press "lift" button again), and engage "spin" button
- sample will sink into the magnet, and spinning will start

### 2. Data acquisition [summary: new → zz → ns → lock → rga → zg ]

- **new <CR>** set-up new file for NMR experiments connected to this sample (nb. more than one experiment can be recording under same file name using different exp. no. extensions)
  - filename: use your lab book ref. plus an appropriate extension letter to indicate origin of sample, e.g. "**PRB502A**"
  - exp. no.: distinguishes between different NMR experiments run on the same sample; it is recommended that you standardize how this index number is used, e.g. 1 for  $^1\text{H}$  NMR, 2 for  $^{13}\text{C}$  NMR, 3 for DEPT135, 4 for COSY
  - proc. no.: 1
  - solvent: choose solvent from drop-down list as appropriate
  - title: use the following format: sample ref., approx. amount, # of scans, solvent, NMR field strength, e.g.  
**"PRB502A 5 mg 16 sc CDC13 300 MHz"**
- **zz <CR>** a macro which brings up a list of common NMR experiments; choose as appropriate (nb. protonstd 14 or 16 refers to desired spectral width (in ppm), use smaller width useless you anticipate unusually downfield signals)
- **ns <CR>** set number of scans as appropriate; for standard  $^1\text{H}$  NMR's and samples of ca 5 mg, 16 scans is ample (nb. make sure ns is a multiple of 4)
- **lock <CR>** automatically sets deuterium lock; choose appropriate solvent from list; machine will take approx. 5 seconds to find and set the lock and indicates when finished

- shimming: while viewing lock display window, adjust "z" and "z<sup>2</sup>" (on axis) to optimize shim set (i.e. make lock line as high as possible); do z first, then z<sup>2</sup>, then go back to z etc. (nb. activate "fine" control as needed); if line goes off screen reduce "lock gain" (can use this function to raise and lower line as desired); when happy with shim, press "std. by" to prevent accidental change of parameters
- **rga <CR>** automatically set receiver gain; machine will pulse sample with rf and monitor strength of NMR signal; procedure takes 5-10 seconds to complete
  - **zg <CR>** "zero go"; command starts data acquisition and will accumulate FID for the desired number of scans  
  
to stop data acquisition early, type "**halt <CR>**", this stops scans *and* saves FID data; *nb. the alternative stop command, "stop <CR>" does not save collected data to the hard disk*

### 3. Data processing

[summary: **ft** → **phase** → **baseline** → **calibrate** → **integrate** → **pick peak**]

- **ft <CR>** to perform Fourier transform on collected FID data; zero-filling can be conducted prior to FT by typing "**si <CR>**" and then doubling number (e.g. from 32 K to 64 K), this operation is optional and improves the cosmetic appearance of peaks making them smoother
- enter phase correction by clicking on icon (or access from processing menu)
  - zero-order: hold mouse pointer over "0" icon and with left mouse button depressed, move mouse up/down to change zero-order phase angle; adjust until line shape of largest peak in spectrum is nicely symmetrical about its base
  - first-order: hold mouse pointer over "1" icon and with left mouse button depressed, move mouse up/down to change first-order phase angle; adjust so that a peak distant from the largest one has good symmetrical shape about its base
- **abs n <CR>** performs automatic baseline correction without auto-integration
- zoom in on desired ref. peak, click the calibrate icon, move calibration line onto the peak and click left mouse button; enter the desired ppm value (see BRG website for values; nb.  $\delta_{\text{H(CHCl}_3\text{)}} = 7.26$  ppm)
- click integration icon; integrate by dragging integrals across desired regions; click save and exit icon when happy with integrals
- pick-peak by dragging box across all signals of interest; click save and exit icon when happy with peak-picks
- **plot <CR>** enter NMR plot software and call-up a BRG standard plot template