

Laboratory #4Week of February 6

Read: Section 9.6, pp. 425-433 of "Optics" by Hecht

Do: 1. Experiment IV.1: Fabry-Perot Interferometer: find Finesse using He-Ne laser  
2. Experiment IV.2: Fabry-Perot Interferometer: Sodium/Mercury doublet


**Experiment IV.1: Fabry-Perot Interferometer: find Finesse using He-Ne laser**

The goal of this lab is to familiarize yourself with the Fabry-Perot Interferometer, to align it properly, and to measure the finesse using the He-Ne laser. The figure below shows a schematic of the Fabry-Perot Interferometer you will use. As with the Michelson Interferometer, you will use a platform for stability. Be sure to use a target when using the Helium-Neon laser. **Do not under any circumstances look into the interferometer output when using the laser** (even though you may when using the lamps).

To align the interferometer, first adjust the two turning mirrors so the laser beam is parallel to the table top (at an approximate height of 13.5 cm) and parallel to but in between the two sets of holes that lead toward the interferometer. This should place the laser beam approximately in the middle of the mirrors. These mirrors are **EXPENSIVE** dielectric coated mirrors, designed for the He-Ne wavelength, with a reflectivity of 95%. Adjust the fixed mirror so that the reflection returns nearly directly back to the laser. There should now be a streak of reflected spots on the target. Adjust the movable mirror so that the streak collapses to a single spot. At this point you may already noticed some interference fringes. To make things more useful, convert the parallel laser input beam into a diverging beam (i.e., a point source) by inserting a diverging lens (-25 mm works well) in front of the interferometer. You should then see a set of concentric circles (a bull's eye pattern). These are the fringes formed by the interference of the multiple beams that are reflected back and forth between the two mirrors of the interferometer. Slight adjustments of the interferometer mirrors will allow you to center this pattern. A small value of  $d$  will provide a larger pattern, which makes the experiment simpler. However, be very careful not to get close enough for the mirrors to touch. These mirrors are quite **expensive**.

The Fabry-Perot interference pattern obeys the same equation as the pattern in the Michelson Interferometer:

$$2d \cos \theta_m = m\lambda,$$



where  $d$  is the spacing of the two mirrors. The new aspect of this interferometer is that multiple beams produce the interference pattern rather than just two as was the case for the Michelson Interferometer. This results in sharper fringes. The Michelson fringes were sinusoidal, whereas the Fabry-

Perot transmission fringes obey the Airy function equation:

$$I(\delta) = \frac{I_0}{1 + \left(\frac{2F}{\pi}\right)^2 \sin^2 \frac{\delta}{2}},$$

where  $F$  is the finesse (not to be confused with the coefficient of finesse) and  $\delta$  is the phase shift, which is given by  $\delta = 2kd\cos\theta$ . The finesse is the ratio of the fringe spacing to the fringe width  $\Delta f$  (FWHM). You can quote the spacing and width in many different ways, e.g., in terms of frequency, wavelength, phase shift, or fringe number. If we work in frequency space, then the separation between adjacent fringes (e.g.  $m$  and  $m+1$ ) is called the Free Spectral Range and is given by  $FSR = c/2d$ , and the finesse is then  $F = FSR/\Delta f$ .

To measure the finesse, scan the movable mirror and record fringes on the oscilloscope. Adjust the mirrors to make the observed fringes as tall and narrow as possible – a very light touch is required. The finesse is a dimensionless number, which you can measure on your oscilloscope as the ratio of the separation of adjacent peaks to the width of the peaks using either seconds or cm as your units. You might observe multiple (2-3) sets of interleaved fringes caused by multiple laser frequencies. If so, turn the laser off/on briefly and watch the fringes move as the laser warms up. This should make it clear which fringes belong to which laser frequency.

### **Experiment IV.2: Fabry-Perot Interferometer: Sodium or Mercury doublet**

In this experiment you will measure the separation of a pair of lines in the sodium or mercury spectrum. The sodium lamp spectrum is dominated by the bright doublet known as the Sodium D-lines at 589.00 and 589.59 nanometers. The mercury lamp emits many different spectral lines, but you will mostly see the yellow lines because they fall within the range of wavelengths that these mirrors reflect (they are designed for Helium-Neon lasers at 633 nm). We also have some filters to select out the yellow lines; their use will make life easier. The mean wavelength of the mercury doublet is 578 nm. The mercury lamp also emits ultraviolet light, so be careful to keep the glass on the front of the lamp (the glass absorbs the UV light) to avoid giving your eyes a sunburn (quite painful) and do not allow any stray light out from the lamp.

The two wavelengths present in the light from the lamp will produce a pair of interference patterns. There will be two sets of concentric circles corresponding to  $\lambda_1$  and  $\lambda_2$ . These wavelengths are only slightly different, which means that patterns will be nearly identical, with only a slight difference in the angular separation between two adjacent fringes. If the two patterns each have a bright fringe at  $\theta = 0$ , then we cannot distinguish one from the other there. As  $\theta$  increases, one set of circles will start to move away from the other until the bright fringes from one wavelength fall exactly half way between the bright fringes of the other wavelength. As  $\theta$  increases further, we will eventually come to another place where the bright fringes overlap. If there are  $N$  fringes of the  $\lambda_1$  pattern between the center and the next location of maximum overlap, then there will be  $N + 1$  fringes of the  $\lambda_2$  pattern (or  $N - 1$  if  $\lambda_1 < \lambda_2$ ). To actually count these numbers we will move the mirror and watch at  $\theta=0$ . In that case, one finds that displacement of the mirror  $\Delta d$  can be written as

$$\Delta d = N \frac{\lambda_1}{2} = (N + 1) \frac{\lambda_2}{2}.$$

Since  $\lambda_1 = \lambda_2 + \Delta\lambda$ , we can derive an equation for the splitting:

$$\Delta\lambda = \frac{\lambda_2}{N}.$$

Place the vapor lamp close to the input of the interferometer. A piece of frosted or ground glass in front of the lamp will help to reduce the intensity and will also provide a uniform background against which to view the fringes. Look into the output port of the interferometer (or use the CCD). You should see a set of concentric fringes as before.

To measure  $\lambda$  (which will be the average of  $\lambda_1$  and  $\lambda_2$ ) you should count 100 fringes. To measure  $\Delta\lambda$ , you should count fringes between overlaps since they are easier to determine than the halfway points. It will be difficult to actually count all the fringes, but you can infer  $N$ , since you know  $\lambda$  and  $\Delta d$ . You should go through several overlaps to improve your precision. If each overlap is hard to determine within 10 fringes and they are separated by 100 fringes, then your error is about 10% (10/100). But if you go through ten overlaps, then your error is only 1% (10/1000). Estimate the error for your experiment.

Equipment needed:

Item	Qty	Source (part #)
Helium-Neon Laser	1	Melles Griot 05 LHP 121
Al mirror	2	Newport 10D10ER.1
Polarizer	2	Edmund A38,396
-25 mm lens	1	Newport KPX043
Sodium vapor lamp	1	
Mercury vapor lamp	1	
Mercury line filter	1	Oriel
3 finger clamp	1	Chem stores
Magnetic base	1	Thor Labs MB175
<u>Fabry-Perot hardware</u>		
Platform for stability	1	Thor Labs (special)
Mounting posts	4	Thor Labs P3
Translation stage	1	Thor Labs MT1-Z8
Mirror mount	2	Thor Labs KM1
Riser block	1	Thor Labs RB2
Base plate	1	Newport BP2
He-Ne mirrors	2	CVI